

**Thermal Waves experiment
as it used to run
(AKA: *where does your data come from?*)**




**YOU NEED TO RUN THIS .PPTX IN 'PRESENTATION MODE'
TO BE ABLE TO WATCH THE VIDEOS**

Thermal Waves – Equipment used in the Lab


HOT BATH (100°C)



COLD BATH (0°C)




PTFE CYLINDER




Thermal Waves – Diagrams Used in the Script


HOT BATH (100°C)



COLD BATH (0°C)



PTFE CYLINDER



This short clip shows the cycling of the cylinder between the hot and cold baths.


- Note that when in the cold bath, the ice-water mixture is constantly stirred to avoid the build up of a warmer layer of water in the vicinity of the PTFE cylinder.
- No stirring is needed in the hot bath due to the rising bubbles from the boiling water.
- The program used for data acquisition produces two beeps:
 - 1) the first one when there's 4s left in the bath
 - 2) the second one when it is time to change the cylinder from one bath to another




Thermal Waves – Data Acquisition Software

These short clips show how the data acquisition software looks like.


(recorded for a square wave period of 2 min. and 20x speed)



Strong temperature transient



Strong temperature transient



Slightly more stable